GAB- U51

Gabapentin Rapid Test Single Panel (Urine)

INTENDED USE

he GAB Rapid Test Single Panel (Urine) is a rapid visual immunoassay for the qualitative, resumptive detection of GAB in human urine specimens at the cut-off concentrations listed below:

arameter	Calibrator	Cut-off (ng/mL)
:AD	Colonastia	2000

INTRODUCTION

iabapentin (GAB) marketed under the brand name Neurontin among others, is a medication used to eat epilepsy, neuropathic pain, hot flashes, and restless leg syndrome. In epilepsy it may be used for ose with partial seizures, It is recommended as one of a number of first line medications for the eatment of neuropathic pain indiabetic neuropathy, post-herpetic neuralgia, and central neuropathic an. The mechanism of the anticonvolusant action of gabapentin has not been fully described. Several ossible mechanisms for pain improvement have been discussed. Though similar in structure to the adogenous neurotransmitted ABA, gabapentin has not been shown to bind to GABA receptors at ancentrations at or below 1 mM. Gabapentin modulates the action of glutamate decarboxylase (GAD) and branched chain aminotransferase (BCAT), two enzymes involved in GABA biosynthesis. In human drat studies, gabapentin was found to increase (ABA biosynthesis, and to increase non-synaptic iABA neurotransmission in vitro Common side effects include sleepiness and dizziness. Serious side flects may include an increased risk of suicide, aggressive behaviour, and drug reaction with semophilia and systemic symptoms. It is unclear if it is safe duringpregnancy or breastfeeding. Lower uses should be used in people with kidney problems. Gabapentin affects the inhibitory eurotransmitter y-aminoburptic acid (GABA).

PRINCIPLE

he GAB Rapid Test Single Panel (Urine) detects Gabapentin through visual interpretation of color evelopment on the Single Panel. Drug conjugates are immobilized on the test region of the membrane uring testing, the specimen reacts with antibodies conjugated to colored particles and precoated on the imple pad. The mixture then migrates through the membrane by capillary action, and interacts with agents on the membrane. If there are insufficient drug molecules in the specimen, the antibody-colored article conjugate will bind to the drug conjugates, forming a colored band at the test region of the iembrane. Therefore, a colored band appears in the test region when the urine is negative for the drug charge mixed the properties of the drug conjugate on the test region for limited antibody binding sites. This will prevent tachment of the antibody-colored particle conjugate to the test region. Therefore, the absence of a slored band at the test region indicates a positive result. The appearance of a colored band at the control gion serves as a procedural control, indicating that the proper volume of specimen has been added and embrane wicking has occurred.

REAGENTS

ach test consists of a reagent Single Panel mounted in a plastic housing. The amount of each antigen id/or antibody coated on the Single Panel is less than 0.001 mg for antigen conjugates and goat uti-rabbit 1gG antibodies, and less than 0.0015 mg for antibody components.

he control zone of each test contains goat anti-rabbit IgG antibody. The test zone of each test contains ug-bovine protein antigen conjugate, and the conjugate pad of each test contains monoclonal anti-drug titlody and rabbit antibody-colored particle complex.

MATERIALS

Materials Provided

Test Panel

Materials Required but Not provided

No contrata de la contrata del contrata de la contrata del contrata de la contrata del contrata de la contrata de la contrata de la contrata del contrata de la contrata del contrata de la contrata del cont

Positive and negative controls

• Timer

Package insert

Centrifuge

PRECAUTIONS

For professional in vitro diagnostic use only

Do not use after the expiration date indicated on the package. Do not use the test if the foil pouch or canister is damaged. Do not reuse tests.

This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore, recommended that these products be treated as potentially infectious, and handled by observing usual safety precautions (e.g., do not ingest or inhale).

Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.

Read the entire procedure carefully prior to testing

- Do not eat, drink or smoke in the area where specimens and kits are handled. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow standard procedures for the proper disposal of specimens. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- · Humidity and temperature can adversely affect results.
- Used testing materials should be discarded in accordance with local regulations.

STORAGE AND STABILITY

- The kit should be stored at 2-30°C until the expiry date printed on the sealed pouch or canister
- The test must remain in the sealed pouch or closed canister until use.
- Do not freeze.
- Kits should be kept out of direct sunlight.
- Care should be taken to protect the components of the kit from contamination. Do not use if there is
 evidence of microbial contamination or precipitation. Biological contamination of dispensing
 equipment, containers or reagents can lead to false results.

SPECIMEN COLLECTION AND STORAGE


- The GAB Rapid Test Single Panel (Urine) is intended for use with human urine specimens only.
- Urine collected at any time of the day may be used.
- · Urine specimens must be collected in clean, dry containers.
- Turbid specimens should be centrifuged, filtered, or allowed to settle and only the clear supernatant should be used for testing.
- Perform testing immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods. Urine specimens may be stored at 2-8°C for up to 2 days. For long term storage, specimens should be kept below -20°C.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Avoid repeated freezing and thawing of specimens.

If specimens are to be shipped, pack them in compliance with all applicable regulations for transportation of etiological agents.

PROCEDURE

Bring tests, specimens, and/or controls to room temperature (15-30°C) before use it if the specimens are refrigerated.

- Remove the test from its sealed pouch and use it as soon as possible. To obtain a best result, the assay should be performed within one hour.
- Hold the Single Panel at the handle with the product name imprints. Do not touch the membrane part of the Single Panel to avoid contamination.
- Dip the test Single Panel vertically in the urine specimen for at least 10-15 seconds. Do not pass
 the maximum line (MAX) on the test Single Panel when immersing the Single Panel. As the test
 begins to work, you will see color move across the membrane.
- 4. Take the Single Panel out of the specimen afterwards and place it on a non-absorbent flat surface. Start the timer and wait for the colored line(s) to appear. The result should be read at 5 minutes. Do not interpret the result after 10 minutes.



INTERPRETATION OF RESULTS

POSITIVE: Only one colored band appears, in the control region (C). No apparent colored band appears in the test region (T).

NEGATIVE: Two colored bands appear on the membrane. One band appears in the control region (C) and another band appears in the test region (T)



NOTE:

The intensity of color in the test region (T) may vary depending on the concentration of analytes present
in the specimen. Therefore, any shade of color in the test region should be considered negative. Note that
this is a qualitative test only, and cannot determine the concentration of analytes in the specimen.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control band failure.

OUALITY CONTROL

- Internal procedural controls are included in the test. A colored band appearing in the control region
 (C) is considered an internal positive procedural control, confirming sufficient specimen volume
 and correct procedural technique.
- External controls are not supplied with this kit. It is recommended that positive and negative
 controls be tested as a good laboratory practice to confirm the test procedure and to verify proper
 test performance.

LIMITATIONS OF THE TEST

- The GAB Rapid Test Single Panel (Urine) is for professional in vitro diagnostic use, and should be only used for the qualitative detection of Gabapentin.
- 2. This assay provides a preliminary analytical test result only. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method by the National Institute on Drug Abuse (NIDA). Clinical consideration and professional judgment should be applied to any test result, particularly when preliminary positive results are indicated.
- There is a possibility that technical or procedural errors as well as other substances and factors may interfere with the test and cause false results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless
 of the analytical method used. Therefore, please preclude the possibility of urine adulteration prior to
 testing.
- 5 A positive result indicates the presence of a Gabapentin only, and does not indicate or measure intoxication.
- A negative result does not at any time rule out the presence of Gabapentin in urine, as they may be present below the minimum detection level of the test.
- 7. This test does not distinguish between Gabapentin and certain medications

PERFORMANCE CHARACTERISTICS

A. Accuracy

The accuracy of the GAB Rapid Test Single Panel (Urine) was compared and checked against commercially available tests with a threshold value at the same cut-off levels. Urine samples taken from volunteers claiming to be non-users were examined under both tests. The results were >95.7% in agreement.

II. Reproducibility

The reproducibility of the GAB Rapid Test Single Panel (Urine) was verified by blind tests performed at four different locations. Samples with GAB concentrations at 50% of the cut-off were all determined to be negative, while samples with GAB concentrations at 200% of the cut-off were all determined to be positive.

C. Precision

Test precision was determined by blind tests with control solutions. Controls with GAB concentrations at 50% of the cut-off yielded negative results, and controls with GAB concentrations at 150% of the cut-off yielded positive results.

D. Specificity

The following tables list the concentrations of compounds (ng/mL) above which the GAB Rapid Test Single Panel (Urine) identified positive results at 5 minutes.

GAB related compounds	Concentration (ng/ml)
Gabapentin	2000
Pregbalin	> 100000

The following compounds yielded negative results up to a concentration of 100 µg/mL:

(-)-Ephedrine Chlorpheniramine Oxalic Acid (+)-Naproxen Creatine Penicillin-G (+/-)-Ephedrine Dextromethorphan Pheniramine 4-Dimethyllaminoantiyrine Dextrorphan tartrate Phenothiazine Acetaminophen Dopamine Procaine Erythromycin Acetone Protonix Albumin Ethanol Pseudoephedrine Amitriptyline Furosemide Quinidine Ampicillin Glucose Ranitidine Aspartame Guaiacol Glyceryl Ether Sertraline Aspirin Hemoglobin Tyramine Benzocaine Ibuprofen Vitamin C (Ascorbic Acid) Bilirubin Imipramine Trimeprazine b-Phenylethyl-amine Isoproterenol Venlafavine Caffeine Lidocaine Methadone Chloroquine

LITERATURE REFERENCES

- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd ed. Davis. Biomedical Publications, 1982.
- Hawks RL, Chiang CN, eds. Urine Testing for Drugs of Abuse. Rockville. Department of Health and Human Services. National Institute on Drug. Abuse; 1986.
- 3. Substance Abuse and Mental Health Services Administration, Mandatory Guidelines for Federal

ce Drug Testing Programs 53 Federal Register, 1988

VI. Drug-analysis technology—pitfalls and problems of drug testing. Clin Chem. 1987 Oct, appl) 33B-40B

AG. Goodman LS, Gilman A, eds. Goodman and Gilman's The Pharmacological Basis of tics. 6th ed. New York. Macmillan, 1980

GLOSSARY OF SYMBOLS

F	Catalog number		Temperature limitation
	Consult instructions for use	LOT	Batch code
]	In vitro diagnostic medical Single Panel	77	Use by
	Manufacturer	2	Do not reuse



Number 1110021770

REVI.0 Effective date: 2017-09-14