

## **Accuracy Check: Alco-Sensor FST**

Using a .038 Alcohol Gas Tank and 1.5 liter-per-minute regulator.





### **Step 1: Attach Regulator**

Follow instructions on the Alcohol Gas Tank to attach the regulator.





Step 2: Purge Regulator by

pressing button/valve on regulator for approx. 8-10 seconds before running your first accuracy check of the day.



**Step 3: Determine Gas Value** by reviewing the elevation chart on the tank to calculate the alcohol concentration for the 0.038 tank at your elevation. The value can also be determined more easily by utilizing a True-Cal device.





**Step 4: Attach mouthpiece**Securely attach disposable FST

mouthpiece to FST instrument.



**Step 5: Access Maintenance Menu** 

Turn on FST by holding down **OFF** button (under thumb) and **ON** button (under index/trigger finger) for a few seconds until screen displays "RCL".







#### **Step 6: Put FST in Accuracy Check**

mode by pressing the ON button (trigger finger) repeatedly until "ACC" is displayed on the screen. Press the OFF button (thumb) once to select this mode.







# **Step 7: Observe Air Blank** The FST displays the temp., then "BLN" for approx. 3 seconds while it performs an automatic Blank Test. (Listen for a <click>.) The display momentarily reads ".000", then flashes "ACC".

Step 8: Attach FST to regulator by attaching the FST mouthpiece to the tube on the tank regulator.



### **Step 9: Perform Accuracy Check**

Press the regulator button for 7 seconds. During the 5<sup>th</sup> second, press & release the FST's **ON** button (trigger finger) to take the gas sample.







#### Step 9 cont'd: Perform Accuracy Check









#### Step 9 cont'd: Additional Notes on Gas Flow

**NOTE 1:** Gas must be flowing through the mouthpiece during the entire 7-second process, therefore, continue pressing the tank regulator button for another 2 seconds AFTER you press & release the **ON** button during the 5<sup>th</sup> second.


**NOTE 2:** If using a 6 liter-per minute regulator, the gas flow will be strong enough to trigger the FST's automatic sampling mechanism, thereby eliminating the need for the Operator to press & release the **ON** button during the 5<sup>th</sup> second.





# Step 10: Observe Reading and Interpret Results

Results within +/- .005 of the elevationadjusted value of the Alcohol Gas Tank are acceptable. The FST is reading accurately.

(see next 3 slides for acceptable readings)









# Step 10 cont'd: Observe Reading and Interpret Results

Results that are **NOT** within +/- .005 of the elevation-adjusted value are out of acceptable tolerance and the Operator must perform a Calibration Adjustment.

(see next 2 slides for unacceptable readings)







# **Step 11: Record Results**

Record every Accuracy Check result in a Calibration Log to document the accuracy of the Alco-Sensor.

#### NOTE:

Results that are **NOT** within +/- .005 of the elevation-adjusted value are out of acceptable tolerance and the Operator must perform a Calibration Adjustment.

For technical assistance and to order calibration supplies:

AlcoPro, Inc. 800 227-9890 www.alcopro.com