



# ToxCup® Reference Guide

The ToxCup® Drug Screen Cup is a one-step test for the rapid qualitative detection of drugs of abuse and adulterants.

**ToxCup® Procedures** -- Refer to package insert for detailed instructions and technical information.

1. **Remove** the test lid from the sealed pouch.



2. **Twist** the test lid securely onto the specimen cup after collection.




3. **Lean** the cup on its side to activate testing.




4. **Adulteration Pads:** After one minute, compare the colors of the adulterant tests with the corresponding color blocks on the color chart included in the kit. **Do not read after 2 minutes as the color may change.**

Color Blocks



5. **Drugs of abuse test:** Read results once the **control band (C)** appears and the membrane clears (in approximately 5 minutes or less). Results are stable for up to an hour after the control bands form. Interpret according to instructions below.

Control Bands



## Interpreting ToxCup® Test Results

### Negative Result

For each test, two colored bands should be observed:

- One in the **CONTROL (C)** region
- One in the **specific TEST** region

The color of the test band may be slightly darker or lighter than the control band. Any band that can be seen visually is a **negative** result.



Example Interpretation:

- AMP: Negative
- COC: Negative
- MET: Negative
- OPI: Negative
- THC: Negative

### Presumptive Positive Result

A colored band at the **CONTROL (C)** region should be observed. When there is no colored band at the specific **TEST** region, the test is presumptive positive for that particular drug.



Example Interpretation:

- AMP: Negative
- COC: **Presumptive Positive**
- MET: Negative
- OPI: Negative
- THC: **Presumptive Positive**

### Invalid Result

When no band appears in the **CONTROL (C)** region, the **TEST** is **invalid** even if there is a band in the test region. Check test procedures, and samples. Repeat the test with a new device.



Example Interpretation:

- AMP: **Invalid**
- COC: Negative
- MET: Negative
- OPI: **Invalid**
- THC: Negative