


iScreen Urine Dip Test

ISCREEN



DRUG TESTING KITS

- Economical dip drug test
- Vivid lines for clear results
- Rapidly screen for a broad range of illicit and prescription drugs
- Multi panel configurations available
- FDA 501(k) cleared

ALCOPRO
DRUG & ALCOHOL TESTING PRODUCTS

1-800-227-9890
www.alcpro.com

DRUG TESTING KITS

iScreen Urine Dip Test



Rapid Testing.

The iScreen Dip Test is a rapid urine screening test that uses monoclonal antibodies to selectively detect elevated levels of specific drugs of abuse. Considered the most economical dip drug test available, the iScreen can be performed without the use of an instrument, stretching your drug testing budget.

Easy-To-Use. Easy-To-Read.

Dip the tip of the test into the urine specimen for ten seconds, then replace the cap and lay flat. Results are ready to be read in as little as 5 minutes. The easy-to-read test results render remarkably prominent lines for negative results, even for notoriously hard-to-read THC panels. These tests also provide a distinct marker line to ensure that you don't dip the tip too deep, flooding the test.

Greater Flexibility.

The iScreen Dip Drug Test offers configurations that can test for a single drug or up to 12 drugs with just one test, giving you greater flexibility in your drug test program. With the iScreen Dip Test, you can choose only the panels required to detect specific drugs of abuse. Limiting the number of panels greatly reduces costs.



Step 1
Remove cover and dip the tip into specimen for 10 sec.



Step 2
Replace cover, lay on a flat surface and read results in 5 min.

ISCREEN DIP DRUG TEST

Item	Description	1 - 24	Box of 25 Tests
8IS01	THC	.99	\$19.75 (.79 ea)
8IS02	Cocaine	.99	\$19.75 (.79 ea)
8IS03	Opiate	.99	\$19.75 (.79 ea)
8IS06	Benzodiazepine	.99	\$19.75 (.79 ea)
8IS21	THC/cocaine	1.89	\$37.25 (1.49 ea)
8IS31	THC/coc/opiate	2.59	\$54.75 (2.19 ea)
8IS42	THC/coc/opi/mAMP	3.39	\$69.75 (2.79 ea)
8IS51	THC/coc/opi/mAMP/pcp	4.39	\$86.25 (3.45 ea)
8IS53	THC/coc/opi/mAMP/amp	4.39	\$86.25 (3.45 ea)
8IS04	Amphetamine	N/A	\$19.75 (.79 ea)
8IS05	Methamphetamine	N/A	\$19.75 (.79 ea)
8IS08	Barbiturate	N/A	\$19.75 (.79 ea)
8IS09	Methadone	N/A	\$19.75 (.79 ea)
8IS10	MDMA (Ecstasy)	N/A	\$19.75 (.79 ea)
8IS11	Oxycodone	N/A	\$19.75 (.79 ea)
8IS12	Buprenorphine	N/A	\$19.75 (.79 ea)
8IS33	THC/coc/mAMP	N/A	\$54.75 (\$2.19 ea)

SPECIFICATIONS

Test Medium: Urine

Format: Dip stick

Cut-off levels:

AMP (amphetamine)	1000 ng/ml
BAR (barbiturates)	300 ng/ml
BZO (benzodiazepines)	300 ng/ml
COC (cocaine)	300 ng/ml
MDMA (Ecstasy)	500 ng/ml
THC (marijuana)	50 ng/ml
MTD (methadone)	300 ng/ml
mAMP (methamphetamine)	1000 ng/ml
OPI (opiates)	300 ng/ml
OXY (oxycodone)	100 ng/ml
PCP (phencyclidine)	25 ng/ml
Bup (buprenorphine)	10 ng/ml

For Information & Orders Call: **1-800-227-9890**
Order 24 Hours A Day At: www.alcopro.com
Fax Orders To: **1-800-655-5834**

ALCOPRO
DRUG & ALCOHOL TESTING PRODUCTS