

TOXCUP

ToxCup



DRUG TESTING KITS

- Easy test procedure: Lean and read
- Test for five drugs
- Results in five minutes
- Minimizes operator exposure to urine
- CLIA-Waived

ALCOPRO
DRUG & ALCOHOL TESTING PRODUCTS

1-800-227-9890
www.alcopro.com

ToxCup



Easy Test Procedure: Lean and Read

The ToxCup uses an ingenious integrated cup design. It builds the test strips into the inner lid of the collection container. The donor provides a urine specimen in the cup, and the operator attaches the leak proof screw-on lid to secure the specimen. When the operator is ready to perform the test, he leans the cup onto two small legs built into the lid. This allows the specimen to contact the test strips and activates the test. The operator reads results in five minutes.

Test for Five Drugs of Abuse Plus Adulterants

The ToxCup comes in several five-panel configurations to test for the most commonly abused drugs of abuse. Each five-panel configuration is also available with four-panel adulteration strips that detect creatinine, pH, oxidants and nitrite.

CLIA-Waived

The ToxCup has CLIA-waived status. This means that ToxCup test results obtained by persons with no laboratory experience and who have only read the testing instructions are comparable to results obtained by trained laboratory technicians. CLIA-waived status allows the ToxCup to be used in a wide variety of settings and reduces the burden otherwise required by CLIA regulations.

SPECIFICATIONS

Test Medium: Urine

Format: Integrated cup with temperature strip

Cut-off levels:

AMP (amphetamine)	1,000 ng/ml
COC (cocaine)	300 ng/ml
OPI (opiates)	300 ng/ml & 2,000 ng/ml
mAMP (methamphetamine)	500 ng/ml & 1,000 ng/ml
PCP (phencyclidine)	25 ng/ml
THC (marijuana)	50 ng/ml

EASE OF USE



Step 1 — Attach leak-proof lid containing test strips.



Step 2 — Lean cup onto built-in legs to activate test.

TOXCUP

Item	Description	I-24	Box of 25
8T51	THC/cocaine/opiates/amp/pcp	\$7.95	\$168.75 (\$6.75 each)
8T51-A	THC/cocaine/opiates/amp/pcp with Adulteration strips	N/A	\$187.50 (\$7.50 each)
8T52	THC/coc/opi/mAMP/pcp	N/A	\$168.75 (\$6.75 each)
8T53	THC/coc/opi/mAMP/amp	N/A	\$168.75 (\$6.75 each)
8T53-A	THC/coc/opi/mAMP/amp with Adulteration strips	N/A	\$168.75 (\$6.75 each) \$187.50 (\$7.50 each)

For Information & Orders Call: **1-800-227-9890**
 Order 24 Hours A Day At: **www.alcopro.com**
 Fax Orders To: **1-800-655-5834**

ALCOPRO
 DRUG & ALCOHOL TESTING PRODUCTS