FASTECTII



DRUG TESTING KITS

- Easy to use dipstick format
- Detects up to four drugs
- Oxycodone and Ecstasy panels available
- Design minimizes operator exposure to urine
- CLIA-Waived

ALCOPRODRUG & ALCOHOL TESTING PRODUCTS

1-800-227-9890 www.alcopro.com

DRUG TESTING KITS

Fastect II



Easy to use Dipstick Format

Dipstick drug test kits are the simplest of all kits to use. Simply dip the tip of the test into the urine specimen for ten seconds and read the test result in five minutes. The Fastect II dipstick test is one of the easiest ways to perform drug testing. The simplicity of the Fastect II dipstick design makes these some of the most economical tests as well.

Detects Up to Four Drugs of Abuse

The Fastect II can test for up to four drugs with one test. The Fastect II is also available in configurations that can test for one, two, and three drugs. This allows customers to choose the only number of panels required to detect the specific drugs their clients are likely to abuse. Limiting the number of panels is the most certain way to lower drug testing costs.

Minimizes Exposure to Urine

The Fastect II design incorporates a sliding cap that minimizes operator exposure to urine. The operator slides the cap up to expose the test strips to perform the test. After wetting the strips in the urine specimen the operator slides the cap to the down position, effectively covering the wetted strips while the test results develop.

CLIA Waived

The Fastect II has CLIA waived status. This means that Fastect II test results obtained by persons with no laboratory experience and who have only read the testing instructions are comparable to results obtained by trained laboratory technicians. CLIA-waived status allows the Fastect II to be used in a wide variety of settings and reduces the burden otherwise required by CLIA regulations.

SPECIFICATIONS	
Test Medium: Urine Format: Dip stick, retractable cap Cut-off levels:	
AMP (amphetamine)	1000 ng/ml
BAR (barbiturates)	300 ng/ml
BZO (benzodiazepines)	300 ng/ml
COC (cocaine)	300 ng/ml
MDMA (Ecstasy)	500 ng/ml
MTD (methadone)	300 ng/ml
OPI (opiates)	300 & 2000 ng/ml
MET (methamphetamine)	500 & 1000 ng/ml
PCP (phencyclidine)	25 ng/ml
THC (marijuana)	50 ng/ml
MOR (morphine)	10 ng/ml

EASE OF USE



Step 1Slide protective cover up to expose tine



Step 2Dip the tips into specimen.



Slide cover down to cover tips, and read test result.

FASTECT II

Item	Description	I - 49	Box of 50 Tests
8F01	THC	\$1.35	\$52.50 (\$1.05 each)
8F02	Cocaine	\$1.35	\$52.50 (\$1.05 each)
8F03	Opiates 300	\$1.35	\$52.50 (\$1.05 each)
8F06	Benzodiazepine	\$1.35	\$52.50 (\$1.05 each)
8F13*	K-2/spice (synthetic THC)	\$6.00	\$250.00 (\$5.00 each)
8F21	THC/cocaine	\$2.75	\$112.50 (\$2.25 each)
8F3 I	THC/coc/opi	\$4.10	\$180.00 (\$3.60 each)
8F42	THC/coc/opi/mAMP	\$5.00	\$199.50 (\$3.99 each)
8F04	Amphetamine		\$52.50 (\$1.05 each)
8F05	Methamphetamine		\$52.50 (\$1.05 each)
8F07	PCP		\$52.50 (\$1.05 each)
8F08	Barbiturate		\$52.50 (\$1.05 each)
8F09	Methadone		\$52.50 (\$1.05 each)
8F10	MDMA (Ecstasy)		\$52.50 (\$1.05 each)
8FII	Oxycodone		\$52.50 (\$1.05 each)
8F12*	Buprenorphine		\$85.00 (\$1.70 each)
8F33	THC/coc/mAMP		\$180.00 (\$3.60 each)
8F43	THC/coc/opi/bzo		\$199.50 (\$3.99 each)

For Information & Orders Call: 1-800-227-9890
Order 24 Hours A Day At: www.alcopro.com
Fax Orders To: 1-800-655-5834