
Total Recall Software

1-800-227-9890

Recal



www.alcopro.com

Total Recall Software

Designed for Corrections and Law Enforcement

The Total Recall software provides a convenient way to upload test results stored in the memory of the Alco-Sensor IV with Memory instruments onto a computer. Total Recall software is designed specifically for agencies that must keep databases of alcohol testing and prepare reports on their alcohol testing program.

Add-on Fields Create Robust Database

The Total Recall software allows the operator to add additional data to each test record after the test results have been uploaded to a PC. By keeping a test log the operator can flesh out the database to create truly meaningful data. The data fields used by the software specifically target the needs of law enforcement and corrections. The software allows the operator to add additional data to each test record for Subject Name, DOB, Drivers License Number, Officer Name, and Case Number.

Uploads test results from Alco-Sensor IV with Memory

The Total Recall software works exclusively with the Alco-Sensor IV with Memory instrument to transfer test results stored in the instrument to a computer. The software includes a cable to connect the Alco-Sensor IV with Memory to the serial port of a PC. Once the Alco-Sensor IV is connected to a PC a few mouse clicks starts the data upload.

Generate Professional-Looking Reports

Total Recall software includes a report generator to print professional-looking pre-formatted reports from the test database. Report filters allow the operator to easily customize reports to meet any request for data. The software applies automatic report formatting so every report comes out presentation-ready.



Graphical Instructions Guide Operator

The Windows based Total Recall program uses intuitive graphical buttons that make the software a snap to use. Most users can start using the software without reading the manual that comes along with the software.

SPECIFICATIONS

Software format: download or CD

Compatible instrument: Alco-Sensor IV with Memory

Minimum computer configuration: Windows 2000, XP; 512 MB Ram, 100 MB hard drive disk space; available serial port (or USB port with approved USB to serial adapter) available printer port.

RELATED PRODUCTS

	Item	Description	Price
THE STATE OF THE S	206-M	Alco-Sensor IV with Memory	\$795.00

Total Recall Software

Item	Description	Price
347	Total Recall Software	\$395.00
	Includes software on CD, Alco-Sensor IV connecting cable, manual.	

For Information & Orders Call: Order 24 Hours A Day At: Fax Orders To:

1-800-227-9890 www.alcopro.com 1-800-655-5834

