

CALIBRATION SUPPLIES

Mini Alcohol Gas Can



ALCOHOL TESTING SUPPLIES

- Use with Alco-Sensor III and Alco-Sensor FST
- Our most economical calibration standard
- Perform Accuracy Checks & Calibration
- Meets needs of halfway houses, shelters, and substance abuse treatment programs

ALCOPRO
DRUG & ALCOHOL TESTING PRODUCTS

1-800-227-9890
www.alcopro.com

Mini Alcohol Can



Use with Alco-Sensor III and Alco-Sensor FST

The Mini Alcohol Cans are ideal for use with the Alco-Sensor III or Alco-Sensor FST when these devices are used as screening devices in substance abuse treatment programs, shelters, and halfway houses. The Mini Alcohol Cans are appropriate for use in any program that does non-evidential alcohol testing. (Agencies that perform evidential quality alcohol testing should use the evidential quality Alcohol Gas Tank.)

Economical Calibration Standard

The Mini Alcohol Can is the most economical non-evidential gas calibration standard available. The Mini Alcohol Cans are ideal for smaller programs that do 10 to 20 Accuracy Checks annually. Typically these programs have only one or two Alco-Sensors, and do Accuracy Checks once a month. The Mini Alcohol Cans provide an affordable solution to performing Accuracy Checks for these low-volume alcohol testing programs.

Perform Accuracy Checks and Calibration Adjustments

Best practice when using any Alco-Sensor is to do periodic (at least monthly) Accuracy Checks to ensure the Alco-Sensor is working correctly. Performing periodic Accuracy Checks gives staff assurance they can rely on the accuracy of Alco-Sensor test results with clients. Most Alco-Sensors hold their accuracy for months or even years. However, when an Accuracy Check shows that the Alco-Sensor is not reading accurately the operator can perform a Calibration Adjustment using the Mini Alcohol Can. Operators can view and download instructions for performing Accuracy Check and Calibration Adjustment procedures at www.alcopro.com.

Specifications

Compatible Instruments: Alco-Sensor III, Alco-Sensor FST

Number of tests per can: 10

Shelf life: Six Months

Evidential Standard: No. Use for screening only

Dimensions: 7" x 2", 2 oz

Gas volume: approximately 2 liters

Alcohol Concentration: Varies from can to can, between 0.080 and 0.120. Each can labeled with specific value

Evidential Quality Gas Tanks

Item	Description	Price
324-04	Alcohol Gas Tank 108 liter (.038)	\$170.00
324-08	Alcohol Gas Tank 108 liter (.082)	\$170.00
324-10	Alcohol Gas Tank 108 liter (.100)	\$170.00
382-04	Alcohol Gas Tank 34 liter (.038)	\$138.50
382-08	Alcohol Gas Tank 34 liter (.082)	\$138.50
370-04	Alcohol Gas Tank 55 liter (.038)	\$149.50
328	Regulator, 1.5 lpm	\$112.00

Mini Alcohol Can

Item	Description	1-4	5-9	10+
303	Mini Alcohol Can	\$37.00	\$35.00	\$33.00

For Information & Orders Call: 1-800-227-9890
 Order 24 Hours A Day At: www.alcopro.com
 Fax Orders To: 1-800-655-5834

ALCOPRO
 DRUG & ALCOHOL TESTING PRODUCTS