




# ALCO-SENSOR FST

Evidential Breath Alcohol Tester

Evidential accuracy in the palm of your hand.

The Alco-Sensor FST is approved by the U.S. Dept. of Transportation as an evidential breath testing device and is trusted worldwide by law enforcement agencies as a preliminary breath tester of choice. Easy-to-use and durable, the Alco-Sensor FST provides reliable results for a more accurate, more defensible testing program.

- Accuracy exceeds +/- .005 at .100 BrAC
- Automatic breath sampling feature captures a deep lung breath sample
- Manual breath sampling is available in shy lung situations or when a subject is uncooperative
- Intuitive testing menus guide the operator through a test sequence, reducing operator errors
- Passive testing mode allows you to screen many people quickly for alcohol and is ideal for large events or roadside screening situations
- One-handed operation is valuable for law enforcement officers



Test subjects every 10 seconds in passive sampling mode.

**GSA** Advantage!®

Available on GSA Advantage for eligible Federal agencies.

**ALCOPRO**  
DRUG & ALCOHOL TESTING PRODUCTS

# SPECIFICATIONS

<b>Sensor Type:</b>	Electrochemical Oxidation	<b>Calibration:</b>	Microprocessor controlled using gas or wet bath simulator
<b>Measurement Range:</b>	.000 - .040 BrAC	<b>Memory:</b>	Last test
<b>Breath Sampling:</b>	Automatic with manual over-ride	<b>Print Capability:</b>	No
<b>Results Response Time:</b>	5 - 10 sec.	<b>Computer Capability:</b>	No
<b>Power:</b>	2 AA batteries (included)	<b>Operating Temperature:</b>	0 - 50° celcius
<b>Test Modes:</b>	Direct (mouthpiece); Passive (sampling cup)	<b>Size:</b>	5" x 2.5" (6 oz)
<b>Warranty:</b>	One year parts and labor	View our device operation, calibration and accuracy check videos online	

## ALCO-SENSOR FST

Item	Standard Product	Intended Markets
207-LE	Alco-Sensor FST, hard-shell case, 25 mouthpieces, battery, sampling cup, drink sniffer, operator manual	Law enforcement, corrections
207-CAN	Alco-Sensor FST, hard-shell case, 25 mouthpieces, battery, sampling cup, drink sniffer, operator manual, mini alcohol gas can	Substance abuse treatment, schools/colleges, non-DOT workplace testing

## RELATED ITEMS AND SUPPLIES

Item	Description
343	Alco-Sensor FST mouthpiece
842	TestSafe™ bacterial viral filter
353	Passive sampling cup
324	Alcohol gas tank, 105L
328-06	Regulator 6 LPM



Detects alcohol in beverages.

## A Few Good Reasons To Trust Us



**Same day shipping**  
Most orders placed before 2:00 pm EST ship same day.



**Knowledgeable Customer Support**  
Chat online or by phone.



**Trusted Supplier more than 35 years.**  
Supporting our customers since 1982.



**Largest inventory**  
We stock more, so you never run out.