



ALCO-SENSOR IV with Memory

Evidential Breath Alcohol Tester

Evidential accuracy in the palm of your hand.

The Alco-Sensor IV with Memory is approved by the U.S. Dept. of Transportation as an evidential breath testing device and is widely used for DOT and non-DOT workplace testing programs. Easy-to-use and durable, the Alco-Sensor IV with Memory provides reliable results for a more accurate, more defensible testing program.

- Evidential accuracy exceeds +/- .005 at .100 BrAC
- Automatic breath sampling feature captures a deep lung breath sample
- Manual breath sampling is available in shy lung situations or when a subject is uncooperative
- Test results display in seconds as a three-digit number
- Intuitive testing menus guide the operator through a test sequence, reducing operator errors
- Stores more than 70 test results in the internal memory
- Automatic :15 wait period countdown following a positive screening test result



ALCOPRO
DRUG & ALCOHOL TESTING PRODUCTS

SPECIFICATIONS

Sensor Type:	Electrochemical Oxidation (fuel cell)	Calibration:	Microprocessor controlled using gas or wet bath simulator
Measurement Range:	.000 - .040 BrAC	Memory:	More than 70 tests
Breath Sampling:	Automatic with Manual override	Print Capability:	Yes (thermal and slip printers)
Results Response Time:	5 - 10 sec. for negative results; 10 - 30 sec. for positive results	Accuracy:	U.S. DOT approved for evidential use
Power:	9-volt battery (included)	Operating Temperature:	0 - 40° celcius
Test Modes:	Direct	Size:	3" x 5" x 1" (11 oz)
Warranty:	Two year parts and labor	Computer:	Yes with software

ALCO-SENSOR IV with Memory

Item	Standard Product
206-M	Hard-shell case, 25 mouthpieces, battery, operator manual, calibration tool
213-SLP	Slip printer, hard-shell case, 100 mouthpieces, printer paper, printer ink ribbon, batteries, printer cables, operator manual, calibration tool
227-P	Thermal printer, hard-shell case, 100 mouthpieces, printer paper, batteries, printer cables, operator manual, calibration tool
227-P-S	Thermal printer, nylon carry case, 100 mouthpieces, printer paper, batteries, printer cables, operator manual, calibration tool, gas tank, regulator

RELATED ITEMS AND SUPPLIES

Item	Description
312	Alco-Sensor IV mouthpiece
842	TestSafe™ bacterial viral filter
382-04	Alcohol gas tank, 34L
328-06	Regulator 1.5 LPM



Item 227-SLP

A Few Good Reasons To Trust Us



Same day shipping
Most orders placed before 2:00 pm EST ship same day.



Knowledgeable Customer Support
Chat online or by phone.



Trusted Supplier more than 35 years.
Supporting our customers since 1982.



Largest inventory
We stock more, so you never run out.