



# ALCO-SENSOR RBT IV

## Evidential Breath Alcohol Tester

Evidential accuracy in the palm of your hand.

The Alco-Sensor RBT IV is approved by the U.S. Dept. of Transportation as an evidential breath testing device and is the most widely used breathalyzer for DOT and non-DOT workplace testing programs. Easy-to-use and durable, the Alco-Sensor RBT IV provides reliable results for a more accurate, more defensible testing program.

- Evidential accuracy exceeds +/- .005 at .100 BrAC
- Automatic breath sampling feature captures a deep lung breath sample
- Manual breath sampling is available in shy lung situations or when a subject is uncooperative
- Automatic :15 wait period countdown following a positive screening test result
- Test results display in seconds as a three-digit number
- Intuitive testing menus guide the operator through a test sequence, reducing operator errors
- Most widely used breathalyzer for DOT testing



Item 211-K

## SPECIFICATIONS

<b>Sensor Type:</b>	Electrochemical Oxidation (fuel cell)	<b>Calibration:</b>	Microprocessor controlled using gas or wet bath simulator
<b>Measurement Range:</b>	.000 - .040 BrAC	<b>Memory:</b>	Last test
<b>Breath Sampling:</b>	Automatic with Manual override	<b>Print Capability:</b>	Yes
<b>Results Response Time:</b>	5 - 10 sec. for negative results; 10 - 30 sec. for positive results	<b>Accuracy:</b>	U.S. DOT approved for evidential use
<b>Power:</b>	9-volt battery (included)	<b>Operating Temperature:</b>	0 - 40° celcius
<b>Test Modes:</b>	Direct	<b>Size:</b>	12" x 9" x 3" (12lbs in case)
<b>Warranty:</b>	Two year parts and labor	<b>RBT IV Printer:</b>	Rechargeable battery good for up to 900 tests

### ALCO-SENSOR RBT IV

Item	Standard Product	Intended Markets
211-K	Alco-Sensor IV, RBT printer, hard-shell case, 100 mouthpieces, printer paper, printer ink ribbon, batteries, operator manual, calibration tool	DOT workplace testing, non-DOT workplace testing
211-K-S	Alco-Sensor IV, RBT printer, nylon case, 100 mouthpieces, printer paper, printer ink ribbon, batteries, operator manual, calibration tool, 34 liter dry gas tank, tank regulator	

### RELATED ITEMS AND SUPPLIES

Item	Description
312	Alco-Sensor IV mouthpiece
842	TestSafe™ bacterial viral filter
382-04	Alcohol gas tank, 34L
328-06	Regulator 1.5 LPM



Item 211-K-S

### A Few Good Reasons To Trust Us



**Same day shipping**  
Most orders placed before 2:00 pm EST ship same day.



**Knowledgeable Customer Support**  
Chat online or by phone.



**Trusted Supplier more than 35 years.**  
Supporting our customers since 1982.



**Largest inventory**  
We stock more, so you never run out.