



# ALCO-SENSOR VXL

## Evidential Breath Alcohol Tester

### Evidential accuracy in the palm of your hand.

The Alco-Sensor VXL is approved by the U.S. Dept. of Transportation as an evidential breath testing device and is widely used for DOT and non-DOT workplace testing. Easy-to-use and durable, the Alco-Sensor VXL provides reliable results for a more accurate, more defensible testing program.

- Evidential accuracy exceeds +/- .005 at .100 BrAC
- Automatic breath sampling feature captures a deep lung breath sample
- Manual breath sampling is available in shy lung situations or when a subject is uncooperative
- Automatic :15 wait period countdown following a positive screening test result
- Test results display in seconds as a three-digit number
- Intuitive testing menus guide the operator through a test sequence, reducing operator errors
- Test subjects without being tethered to a printer for optimal mobility
- Stores thousands of test results in the internal memory



**ALCOPRO**  
DRUG & ALCOHOL TESTING PRODUCTS

## SPECIFICATIONS

<b>Sensor Type:</b>	Electrochemical Oxidation (fuel cell)	<b>Calibration:</b>	Microprocessor controlled using gas or wet bath simulator
<b>Measurement Range:</b>	.000 - .040 BrAC	<b>Memory:</b>	More than 1,000 tests
<b>Breath Sampling:</b>	Automatic with Manual override	<b>Print Capability:</b>	Yes (with docking station and thermal printer)
<b>Results Response Time:</b>	3 - 5 sec. for negative results; 5 - 10 sec. for positive results	<b>Accuracy:</b>	U.S. DOT approved for evidential use
<b>Power:</b>	4 AA batteries (included)	<b>Operating Temperature:</b>	-5° to 50° celcius
<b>Test Modes:</b>	Direct	<b>Size:</b>	7" x 3" x 2 1/8" (12 oz)
<b>Warranty:</b>	Two year parts and labor	<b>Computer:</b>	Yes with docking station and Accelerator software

### ALCO-SENSOR VXL

Item	Standard Product	Intended Markets
208-DOT	Alco-Sensor VXL, hard-shell case, 25 mouthpieces, batteries, operator manual	DOT workplace testing, Non-DOT workplace testing
225	Alco-Sensor VXL, slip printer, hard-shell case, 100 mouthpieces, batteries, printer paper, printer ink ribbon, docking station, manual	

### RELATED ITEMS AND SUPPLIES

Item	Description
380	Alco-Sensor VXL mouthpiece
842	TestSafe™ bacterial viral filter
382-04	Alcohol gas tank, 34L
328-06	Regulator 1.5 LPM



Item 225

### A Few Good Reasons To Trust Us



**Same day shipping**  
Most orders placed before 2:00 pm EST ship same day.



**Knowledgeable Customer Support**  
Chat online or by phone.



**Trusted Supplier more than 35 years.**  
Supporting our customers since 1982.



**Largest inventory**  
We stock more, so you never run out.