



iCup

Rapid Drug Test Kit: Integrated Cup

Easy to Use. Rapid Results.

The iCup has a drug detection card built into the urine specimen container, minimizing contact with the urine specimen. No need to pipette the urine onto a test card. Results are ready in 5 minutes.

The unique flat front design enables you to photocopy the test results from the test card along with a provided report template, to create a permanent test record with donor information.

Send the iCup container with the specimen to a laboratory to confirm non-negative results.

- Multiple drug panel configurations are available
- Easy to read results, just peel back the label
- Rapid results in 5 minutes
- Temperature strip helps to validate the specimen
- Flat front panel allows you to photocopy test results
- Secure leak-proof lid keeps the urine securely in the cup
- Aldulteration panel available on 4, 5, 8 and 10-panel test kits



iCup RX screens for prescription drugs of abuse.

SPECIFICATIONS

Test Medium:	Urine	Format:	Integrated cup with temperature stripe
Cut-Off Levels			
AMP (Amphetamines)	1,000 ng/mL	BAR (Barbiturates)	300 ng/mL
BZO (Benzodiazepines)	300 ng/mL	COC (Cocaine)	300 ng/mL
MDMA (Ecstasy)	500 ng/mL	THC (Marijuana)	50 ng/mL
MTD (Methadone)	300 ng/mL	mAMP (Methamphetamines)	1,000 ng/mL
OPI (Opiates)	2,000 ng/mL	OXY (Oxycodone)	100 ng/mL
PCP (Phencyclidine)	25 ng/mL	TCA (Tricyclic Antidepressants)	1,000 ng/mL

iCup and iCup RX (*CLIA Waived)

Item	Drug Testing Panels
8i41*	4-Panel Test (THC/COC/mAMP/OPI)
8i53*	5-Panel Test (THC/COC/AMP/mAMP/OPI)
8i81	8-Panel Test (THC/COC/AMP/mAMP/OPI/PCP/BZO/BAR)
8i101	10-Panel Test (THC/COC/AMP/mAMP/OPI/PCP/BZO/BAR/MTD/TCA)
8iRX	iCup RX 5-Panel Test (BUP/BZO/MTD/OPI/OXY)

EASE OF USE



Step 1
Collect the specimen and check the temperature.



Step 2
After 5 minutes, peel back the label to read results.



Step 3
The leak-proof lid allows you to photocopy test results.

A Few Good Reasons To Trust Us



Same day shipping
Most orders placed before 2:00 pm EST ship same day.



Knowledgeable Customer Support
Chat online or by phone.



Trusted Supplier more than 35 years.
Supporting our customers since 1982.



Largest inventory
We stock more, so you never run out.