

The sweat patch has a successful history against court challenges. High Times magazine gave an unintentional endorsement when it advised readers to avoid wearing the patch, if possible.



# Sweat Patch

## Drug Test Kits


### Continuous and Reliable Drug Monitoring.

The Sweat Patch is designed to be worn for seven days and is available as a standard six-panel or with extended opiate panel.

Each pad consists of a collector pad and strong adhesive film. If the subject uses drugs while wearing the patch, a small portion of the drug and/or drug metabolite is excreted into the sweat and absorbed by the collector pad.

Common adulteration methods are ineffective with the sweat patch. Each patch is analyzed by a certified laboratory with non-negative results confirmed using the GC/MS methodology to produce evidential-quality results.

- Resistant to tampering and adulteration attempts
- Includes laboratory screening and GS/MS confirmation of non-negatives
- For court-ordered treatment and drug diversion programs
- Continuous monitoring up to 7 days
- Water resistant so it stays on during showers and swimming



# SPECIFICATIONS

**Test Medium:** Sweat **Format:** Absorbent collector pad

**Cutoff Levels:** (screening cutoff levels are the same to cutoff levels in urine testing)

THC: 1.5 ng/ml	PCP: 7.5 ng/ml
Cocaine: 10 ng/ml	Hydrocodone: 10 ng/ml
Opiates: 10 ng/ml	Hydromorphone: 10 ng/ml
Amphetamine: 10 ng/ml	Oxycodone: 10 ng/ml
Methamphetamine: 10 ng/ml	Oxymorphone: 10 ng/ml

## SWEAT PATCH DRUG TEST

Item	Product
882	Standard Drug Test Panel (amphetamines, cocaine, hydrocodone, hydromorphone, methamphetamines, opiates, PCP, THC)
882-FT	Standard Drug Test Panel Plus Fentanyl Panel (amphetamines, cocaine, hydrocodone, hydromorphone, methamphetamines, opiates, PCP, THC, fentanyl)
882-Opiate	Extended Opiate Panel (amphetamines, benzodiazepine, buprenorphine, cocaine, fentanyl, hydrocodone, hydromorphone, methadone, methamphetamines, opiates, oxycodone, oxymorphone, PCP, THC)
882-Overlay	Sweat Patch Overlay (protective overlay recommended for subjects who perform rigorous work or workouts)

## RELATED ITEMS AND SUPPLIES

Item	Description
860	Nitrile gloves (sm, med, lg, xl)
861	Latex gloves (sm, med, lg, xl)



Overlay protects patch

## A Few Good Reasons To Trust Us



**Same day shipping**  
Most orders placed before 2:00 pm EST ship same day.



**Knowledgeable Customer Support**  
Chat online or by phone.



**Trusted Supplier more than 40 years.**  
Supporting our customers since 1982.



**Largest inventory**  
We stock more, so you never run out.