



# iCup

## Rapid Drug Test Kit: Integrated Cup


### Easy to Use. Rapid Results.

The iCup has a drug detection card built into the urine specimen container, minimizing contact with the urine specimen. No need to pipette the urine onto a test card. Results are ready in 5 minutes.

The unique flat front design enables you to photocopy the test results from the test card along with a provided report template, to create a permanent test record with donor information.

Send the iCup container with the specimen to a laboratory to confirm non-negative results.

- Multiple drug panel configurations are available
- Easy to read results, just peel back the label
- Rapid results in 5 minutes
- Temperature strip helps to validate the specimen
- Flat front panel allows you to photocopy test results
- Secure leak-proof lid keeps the urine securely in the cup
- Aldulteration panel available on 4, 5, 8 and 10-panel test kits



## SPECIFICATIONS

<b>Test Medium:</b>	Urine	<b>Format:</b>	Integrated cup with temperature stripe
AMP (Amphetamines)	1,000 ng/mL	BAR (Barbiturates)	300 ng/mL
BZO (Benzodiazepines)	300 ng/mL	COC (Cocaine)	300 ng/mL
MDMA (Ecstasy)	500 ng/mL	THC (Marijuana)	50 ng/mL
MTD (Methadone)	300 ng/mL	mAMP (Methamphetamines)	1,000 ng/mL
OPI (Opiates)	2,000 ng/mL	OXY (Oxycodone)	100 ng/mL
PCP (Phencyclidine)	25 ng/mL	TCA (Tricyclic Antidepressants)	1,000 ng/mL

Item	Drug Testing Panels with Adulteration Panels (*CLIA Waived)
8i41*	4-Panel Test (THC/COC/mAMP/OPI)
8i53/8i51*	5-Panel Tests (THC/COC/AMP/mAMP/OPI) AND (THC/AMP/COC/OPI/PCP)
8i61	6-Panel Test (THC/AMP/BZO/COC/mAMP/OPI)
8i81	8-Panel Test (THC/COC/AMP/mAMP/OPI/PCP/BZO/BAR)
8i101	10-Panel Test (THC/COC/AMP/mAMP/OPI/PCP/BZO/BAR/MTD/TCA)
8i111	11-Panel Test (THC/AMP/BAR/BZO/COC/mAMP/MTD/OPI/OXY/PCP/PPX)
8i112	12-Panel Test (THC/AMP/BAR/BZO/COC/mAMP/MTD/OPI/OXY/PCP/PPX/TCA)

## EASE OF USE



### Step 1

Collect the specimen and check the temperature.



### Step 2

After 5 minutes, peel back the label to read results.



### Step 3

The leak-proof lid allows you to photocopy test results.

## A Few Good Reasons To Trust Us



### Same day shipping

Most orders placed before 2:00 pm EST ship same day.



### Knowledgeable Customer Support

Chat online or by phone.



### Trusted Supplier more than 40 years.

Supporting our customers since 1982.



### Largest inventory

We stock more, so you never run out.