




HOUND® CANNABIS BREATHALYZER

COLLECT + SEND

The HOUND® CANNABIS BREATHALYZER – Collect + Send is the first breath test designed for the era of cannabis legalization – limiting detection to the workday. This novel technology targets only THC molecules, automates breath collection, and enables ultra-sensitive lab-based analysis.

WHY CHOOSE THE HOUND® CANNABIS BREATHALYZER?

- + On-site collection with lab-based processing of results
- + Standardizes and automates breath collections to optimize validity
- + Single-use, tamper-evident cartridges yield no-mess collections
- + Smart base station with lab integration facilitates paperless management and electronic chain of custody
- + Integration with the Hound Labs Retriever™ results portal provides convenient solution to order, track, and manage test results in real-time



WHY TEST BREATH FOR CANNABIS?

- + Renders shortest cannabis detection window – within a few hours of the test
- + Legacy methods detect cannabis use for days, weeks, or months after use
- + Supports policies that deter cannabis use during the workday
- + Promotes safety, protects employee privacy, encourages equity
- + Reduces risk and improves employee recruitment and retention
- + Aligns with emerging cannabis laws such as California A.B. 2188 and Washington S.B. 5123 which limit the use of traditional cannabis tests that extend beyond the workday

ABOUT HOUND LABS

Hound Labs, Inc. is committed to developing ultra-sensitive, automated, and analytical breath technologies that will pave the way to practical and fair solutions for leading public health and safety issues.

The HOUND® CANNABIS BREATHALYZER is intended to detect recent cannabis use. It does not measure whether, or how much, a person is impaired. It is intended solely for use in employment, law enforcement, and insurance settings. It should not be used for any medical or therapeutic purposes, or for any Federal drug testing programs, such as programs run by the Substance Abuse and Mental Health Services Administration (SAMHSA), the Department of Transportation (DOT), and the U.S. military.

SAFETY + FAIRNESS™

July 2023

Call AlcoPro at **800.227.9890** for a product demonstration or quote.



HOUND® CANNABIS BREATHALYZER

COLLECT + SEND

HOW DOES IT WORK?

1 ATTACH

Standardized samples help prevent invalid results


2 BREATHE

Automated breath collections reduce human error

3 PROCESS

Collect + Send results

The HOUND® CANNABIS BREATHALYZER – Collect + Send is designed to be intuitive and easy to use for employers, collectors, and donors. A smart base station provides simple touchscreen instructions to guide users through the collection process.



ADDITIONAL SPECIFICATIONS

Base station: 5.5"H x 14.25"W x 9.25"D / ~10.4 lbs

Cartridge: 4.5"H x 1.5"W (compressed 3.25"H) / ~2.5oz

The HOUND® CANNABIS BREATHALYZER is intended to detect recent cannabis use. It does not measure whether, or how much, a person is impaired. It is intended solely for use in employment, law enforcement, and insurance settings. It should not be used for any medical or therapeutic purposes, or for any Federal drug testing programs, such as programs run by the Substance Abuse and Mental Health Services Administration (SAMHSA), the Department of Transportation (DOT), and the U.S. military.

SAFETY + FAIRNESS™

July 2023

ALCOPRO.COM/HOUND-CANNABIS-BREATHALYZER