

Choosing an Evidential Breath Alcohol Instrument



Do you plan to take disciplinary action for positive test results AND/OR need an instrument with legally defensible test results?

if no

if yes

Non-Evidential Alcohol Testing: Screening Tests using the Alco-Screen or Alco-Screen 02 saliva screening tests

Evidential Breath Alcohol Testing: Do you conduct DOT Tests OR Non-DOT Tests?...OR does your company policy mimic DOT testing regulations?

if non-dot

*if dot or mimics dot

if no

if yes

Do you need to print screening test results?

Do you need a portable instrument capable of operating without AC power?

if no

if yes

DOT Versions

- Alco-Sensor VXL with Slip Printer

DOT Versions

- Alco-Sensor RBT VXL (Data Lite and Full Data systems)

Generic/Non-DOT Versions

- Alco-Sensor FST

Generic/Non-DOT Versions

- Alco-Sensor VXL Universal Screener

*DOT testing is performed when you have employees performing safety-sensitive duties which are covered under the jurisdiction of the U.S. Dept. of Transportation and using Federal Regulations 49 CFR Part 40.

Training Considerations:

- If performing DOT testing...DOT Breath Alcohol Technician Training is required
- If performing Non-DOT testing...Breath Alcohol Technician Operator Training is recommended

www.alcopro.com

e: info@alcopro.com

p: 800.227.9890

f: 800.655.5834

When you're responsible for drug and alcohol testing, turn to the company with more than 40 years of experience delivering quality products, expert training, and fast, reliable service.

