



ALCO-SENSOR VXL

Evidential Breath Alcohol Tester

Evidential accuracy in the palm of your hand.

The Alco-Sensor VXL is approved by the U.S. Dept. of Transportation as an evidential breath testing device and is widely used for substance abuse treatment, probations, and non-DOT workplace testing. Easy-to-use and durable, the Alco-Sensor VXL provides reliable results for a more accurate, more defensible testing program.

- Evidential accuracy exceeds +/- .005 at .100 BrAC
- Automatic breath sampling feature captures a deep lung breath sample
- Manual breath sampling is available in shy lung situations or when a subject is uncooperative
- Test results display in seconds as a three-digit number
- Intuitive testing menus guide the operator through a test sequence, reducing operator errors
- Passive test mode allows the operator to test multiple subjects quickly
- Test subjects without being tethered to a printer for optimal mobility
- Stores thousands of test results in the internal memory



Item 208

ALCOPRO
DRUG & ALCOHOL TESTING PRODUCTS

SPECIFICATIONS

Sensor Type:	Electrochemical Oxidation (fuel cell)	Calibration:	Microprocessor controlled using gas or wet bath simulator
Measurement Range:	.000 - .040 BrAC	Memory:	More than 1,000 tests
Breath Sampling:	Automatic with Manual override	Print Capability:	Yes (with docking station and thermal printer)
Results Response Time:	3 - 5 sec. for negative results; 5 - 10 sec. for positive results	Accuracy:	U.S. DOT approved for evidential use
Power:	4 AA batteries (included)	Operating Temperature:	-5° - 50° celcius
Test Modes:	Direct and Passive	Size:	7" x 3" x 2 1/8" (12 oz)
Warranty:	One year parts and labor	Computer:	Yes with docking station and Accelerator software

ALCO-SENSOR VXL

Item	Standard Product	Intended Markets
208	Hard-shell case, 25 mouthpieces, batteries, drink sniffer, 2 passive sampling cups, operator manual	Substance Abuse Treatment, Law Enforcement, Non-DOT workplace testing
208-BAR	Built-in barometer, hard-shell case, 25 mouthpieces, batteries, drink sniffer, 2 passive sampling cups, operator manual	
224-KIT	Thermal printer, nylon carry case, 100 mouthpieces, batteries, docking station, printer paper, drink sniffer, 2 passive sampling cups, manual	

RELATED ITEMS AND SUPPLIES

Item	Description
380	Alco-Sensor VXL mouthpiece
842	TestSafe™ bacterial viral filter
382-08	Alcohol gas tank, 34L
328-06	Regulator 1.5 LPM



Item 224-KIT

A Few Good Reasons To Trust Us



Same day shipping
Most orders placed before 2:00 pm EST ship same day.



Knowledgeable Customer Support
Chat online or by phone.



Trusted Supplier more than 35 years.
Supporting our customers since 1982.



Largest inventory
We stock more, so you never run out.