



iScreen Nicotine Test

Instant Cotinine Oral Fluid Screening Test

Detect Recent Nicotine Use in Minutes

The iScreen® Oral Fluid Nicotine Test provides a fast, reliable way to screen for recent nicotine and tobacco use without the need for urine collection or invasive procedures. Designed for workplace wellness programs, healthcare providers, insurance screening, research studies, and tobacco cessation programs, this easy-to-use test delivers clear pass/fail results in just 10 minutes.

Using a built-in saliva collection swab, the all-in-one device simplifies specimen collection while reducing handling steps and improving user convenience.

Product Features

- Integrated saliva collection swab
- Clear pass/fail visual results in 10 minutes
- Approximately one to two day detection window
- Detects cotinine at a 30 ng/mL cutoff
- Tamper-evident laboratory confirmation capability
- Simple three-step testing process



ALCOPRO
DRUG & ALCOHOL TESTING PRODUCTS

Why Choose the iScreen Oral Fluid Nicotine Test?

Fast Results

Receive easy-to-read positive or negative nicotine screening results in 10 min.

Non-Invasive Sample Collection

The integrated oral fluid collection swab eliminates the need for urine collection, making the test process easier to administer.

Detects Recent Nicotine Use

The test qualitatively detects cotinine - the primary metabolite of nicotine - at 30 ng/mL cutoff, identifying nicotine use typically within the previous 1-2 days.

Laboratory Confirmation Ready

When confirmation testing is required, simply apply a tamper-evident seal and send the device directly to the laboratory.

Built-In Quality Control

An integrated procedural control verifies the test has functioned correctly and results can be interpreted with confidence.

Item	Standard Product	Intended Markets
8S019	Instant oral fluid nicotine test kit; sold individually or in boxes of 25 test kits (quantity pricing is available)	Workplace testing insurance rate-reduction programs; Health or smoking cessation programs



1. Insert sponge end into the mouth and swab the inside of the mouth and tongue for at least three minutes.
2. Insert the collector vertically into the cap and press down firmly until the cap reaches the thread. Twist the handle clockwise 180 degrees to tighten the cap until the red arrow lines up with the other white mark.
3. Place the device horizontally on a flat surface (result side up), and read results at 10 minutes.

A Few Good Reasons To Trust Us



Same day shipping

Most orders placed before 2:00 pm EST ship same day.



Knowledgeable Customer Support

Chat online or by phone.



Trusted Supplier more than 40 years.

Supporting our customers since 1982.



Largest inventory

We stock more, so you never run out.